


TYPE EVM-D1.2



Pumpdown time for 20 ltr vessel



Technical Data		TYPE EVM-D1.2
		Twin head vacuum
Nominal displacement	m ³ /h	1,2
Ultimate pressure	mbar	50
Motor power	W	250
Material Diaphragm & Valves		HNBR
Material Pump Head		PTFE coated
Weight	kg	10
Noise level	dB	< 50
Admissible ambient temp.	°C	10 to 40
Connections	mm	6
Dimensions LxWxH	mm	282x130x210

Ordering Information	TYPE EVM-D1.2
Single phase motor 220-240V, 50 Hz (IEC)	623.000
Accessories	See page 8