


TYPE EVM-S0.7



Pumpdown time for 5 ltr vessel



| Technical Data | | TYPE EVM-S0.7 | |
|---|-------------------|--------------------|--|
| | | Single head vacuum | |
| Nominal displacement | m ³ /h | 0,72 | |
| Ultimate pressure | mbar | 300 | |
| Motor power | W | 75 | |
| Pump head | | Nylon | |
| Weight | kg | 4,0 | |
| Noise level | dB | < 50 | |
| Admissible ambient temp. | °C | 10 to 40 | |
| Connections | mm | 6 | |
| Dimensions LxWxH | mm | 195x98x156 | |
| Ordering Information | | TYPE EVM-S0.7 | |
| Single phase motor 220-240V, 50 Hz (IEC) | | 621.000 | |
| Accessories | | See page 8 | |