


TYPE EVM-S1.2



Pumpdown time for 10 ltr vessel



Technical Data		TYPE EVM-S1.2	
		Single head vacuum	
Nominal displacement	m ³ /h	1,2	
Ultimate pressure	mbar	250	
Motor power	W	200	
Material Diaphragm & Valves		HNBR	
Material Pump Head		PTFE coated	
Weight	kg	7,5	
Noise level	dB	< 50	
Admissible ambient temp.	°C	10 to 40	
Connections	mm	6	
Dimensions LxWxH	mm	235x140x210	

Ordering Information		TYPE EVM-S1.2	
Single phase motor 220-240V, 50 Hz (IEC)		622.000	
Accessories		See page 8	