


TYPE EVM-S1.8



Pumpdown time for 10 ltr vessel



Technical Data		TYPE EVM-S1.8
		Single head vacuum/pressure
Nominal displacement	m ³ /h	1,8
Ultimate pressure	mbar	250
Motor power	W	200
Material Diaphragm & Valves		HNBR
Material Pump Head		PTFE coated
Weight	kg	8,0
Noise level	dB	< 50
Admissible ambient temp.	°C	10 to 40
Connections	mm	6
Dimensions LxWxH	mm	215x165x270

Ordering Information	TYPE EVM-S1.8
Single phase motor 220-240V, 50 Hz (IEC)	624.000
Accessories	See page 8