

Accessories for EVM series



Syringe Filters

Disposable syringe filters designed to provide fast and efficient filtration of aqueous and organic solutions. They are made of a wide variety of different membranes filters using the most advanced methods and design features available today. These syringe filters are ideal for numerous applications in pharmaceutical, environmental, biotechnology, food/beverage and agricultural testing laboratories. The syringe filters are composed of a pure polypropylene housing, and are heat sealed without the use of glues or sealants. On your demand we can produce the syringe filters in sterilized package and ultrasonic welded.

Applications:

- Protein and enzyme filtration, sterilization
- Biological fluid filtration sterilization
- Tissue culture media sterilization
- Diagnostic cytology
- Receptor binding studies
- Enhanced recovery of fastidious gram positive organisms

Ordering Information			Syringe Filter Set 100 pieces per set
Ø25 mm	PTFE	0,45 µm	629.950
Ø25 x 0,2mm	PES	0,45 µm	629.951
Ø25 x 0,2mm	Nylon6	0,45 µm	629.952
Ø25 x 0,2mm	Nylon66	0,45 µm	629.953
Ø25 x 0,2mm	PVDF	0,45 µm	629.954
Ø13 mm	PTFE	0,45 µm	629.955
Ø13 x 0,2mm	PES	0,45 µm	629.956
Ø13 x 0,2mm	Nylon66	0,45 µm	629.957
Ø13 x 0,2mm	PVDF	0,45 µm	629.958